## pH Measurements/Indicator Kit #81R Student Work Sheet and Guide

An important method of determining pH values involves the use of "indicators". Indicators are certain organic substances that have the property of changing color in dilute solutions when the hydrogen-ion concentration of the solution reaches a definite value. For example, phenolphthalein solution is a colorless liquid in any aqueous solution of which the pH is less than 9. In solutions for which the pH is greater than 9, the phenolphthalein imparts a red or pink color to the solution. Substances like phenolphthalein are called acid-base indicators and are often used for determining the approximate pH of solutions.

This lab activity is designed to allow you to study and gain a better understanding of pH and its relationship to acids and bases as well as its effect on indicator dyes.

You will be working individually or in a team depending on the instructions given to you by your teacher.

One CHEMPLATE® per student/team is required. Several paper towels are needed.

In all activities replace the cap after use and do not touch the tip of the bottle to the CHEMPLATE \* or to the solution in the cavities. When dropping, hold the bottle vertically then squeeze gently.

## **ACTIVITY I. THE EFFECT OF pH ON INDICATOR DYES**

**PURPOSE:** To determine the color change for each of three indicators: methyl orange, brom thymol blue and phenolphthalein. **PROCEDURE:** 

1. Place the CHEMPLATE® on a sheet of white paper.

1. Methyl orange solution is what color?

- 2. Place one drop of methyl orange in cavities #1 and #2.
- 3. Place one drop of brom thymol blue in cavities #5 and #6.
- **4.** Place one drop of phenolphthalein in cavities #9 and #10.
- 5. Carefully add one drop of pH 1 to cavities #1, #5 and #9.
- 6. Carefully add one drop of pH 13 to cavities #2, #6 and #10.

## **OBSERVATIONS:**

2. Methyl orange changes to		in an acid (pH 1) solution.
3. Methyl orange changes to	(color)	in a basic (pH 13) solution.
4. Brom thymol blue solution is what color?	(color)	
5. Brom thymol blue changes to		in an acid (pH 1) solution.
6. Brom thymol blue changes to	(color)	in a base (pH 13) solution.
7. Phenolphthalein solution is what color?	(color)	
8. Phenolphthalein changes to		in an acid (pH 1) solution.
9. Phenolphthalein changes to	(color)	in a base (pH 13) solution.
Rinse the CHEMPLATE® with tap water and dr		
ACTIVITY II. DETER PURPOSE: To establish the specific pH range PROCEDURE: IIA.	and the second s	NGE with INDICATOR DYES unges for each indicator.
1. Place 1 drop of methyl orange in each	cavity #'s 1-7.	
		cavity #2; 1 drop of pH 5 to cavity #3; 1 drop of pH 7 ocavity #6; And 1 drop of pH 13 to cavity #7.

11B.	
3. Place 1 drop of brom thymol blue in each cavity #'s 1-7	
4. Repeat step 2 from Procedure IIA.	
OBSERVATIONS:	
brom thymol blue changes to	atpH.
Rinse the CHEMPLATE <sup>®</sup> and dry. (color)	(number)
IIC.	
5. Place 1 drop of phenolphthalein in each cavity #'s 1-7.	
6. Repeat step 2 from Procedure IIA.	
OBSERVATIONS:	
Phenolphthalein changes to	at pH.
Rinse the CHEMPLATE® and dry. (color)	(number)
ACTIVITY III. USING UN	NIVERSAL INDICATOR
PROCEDURE:	
1. Place 1 drop of Universal Indicator in each cavity #'s 1-	<b>-7.</b>
2. Repeat step 2 from Procedure IIA.	
OBSERVATIONS:	
1. Describe the color in each cavity.  #1 (pH 1)  #2 (pH 3)  #3 (pH 5)  #4 (pH 7)  #5 (pH 9)  #6 (pH 11)  #7 (pH 13)	
DO NOT MIX OR RINSE THE SOLUTIONS IN THE CHEMPLAT	<b>TE</b> ®
ACTIVITY IV. DETERMINING THE PROCEDURE:	pH OF UNKNOWN SOLUTIONS
1. Place 1 drop of Universal Indicator in cavities #10, #1	l1 and #12.
2. Place 1 drop of Unknown Lin cavity #10.	
3. Place 1 drop of Unknown II in cavity #11.	
4. Place 1 drop of Unknown III in cavity #12.	
5. Compare the colors in cavities #10, #11, and #12 with	the colors in cavities #ls 1.7
OBSERVATIONS:	αιο σοιοίο το σαντίτεο π ο 177.
1. The pH of Unknown I is	
2. The pH of Unknown II is	
3. The pH of Unknown III is	
Student's Name	